

University of Stuttgart

Institute of Combustion and Power Plant Technology

Prof. Dr. techn. G. Scheffknecht

Tar and sulfur measurement and monitoring for biogenic residue gasification

IEA Task 33 Workshop
Gas cleaning, experiences, new developments,
analytics and diagnostics
6th June 2019, Karlsruhe

M. Schmid,
G. Scheffknecht

Agenda



- gasification processes at IFK
- FID online tar analyzer TTA300
- sulfur measurement at IFK
- liquid tar sampling / tar protocol at IFK
- planned new online GC for tar and impurity measurement
- summary and conclusion

Motivation

- mitigation of climate change: to abandon the utilization of fossil carbon




- on contrast: continuing demand of carbon-based products
 - plastics, chemicals
 - fuels for transport (eg. kerosene), fuels for flexible heat and power production

→ utilization of carbon from biomass AND recovery of carbon from waste materials




gasification as conversion process:



- Fluidized bed gasification → high efficiency, high feedstock flexibility
 - current research at IFK:
 - sorption enhanced gasification (SEG)
 - steam-oxygen gasification

Fluidized bed Gasification Processes at IFK


SEG



(Oxy-)SEG (Sorption Enhanced Gasification)

- in-situ CO_2 capture: Limestone shifts CO_2 vom gasifier to regenerator
- H_2 rich and CO_2 lean syngas
- CO_2 -rich off-gas for storage


Steam-oxygen gasification



- autothermal single bed gasification
- syngas suitable for synthesis
- applicable for high-ash containing fuels such as dried sewage sludge


Fluidized bed research and pilot facilities (5-200 kW_{th})

1.5 kW_{th} (el. heated)



Height	0.65 m
Fuel rate	Up to 1 kg/h
Fluidization	BFB

20 kW_{th} (el. heated)



Height	3,5-12 m
Fuel rate	1-5 kg/h
Fluidization	BFB/CFB

200 kW_{th} (not el. heated)




Height	6-10 m
Fuel rate	20-50 kg/h
Fluidization	BFB/CFB

FID Online Total Tar Analyzer TTA300

FID Online Total Tar Analyzer TTA300

produced by Ratfisch Analysensysteme





- measurement total tar concentration
- semi-continuous measurement every 1,5 min
- device operation temperature 300°C
- state of the art software and easy operation
- suitable for tar concentration $> 1 \text{ g/m}^3$
→ raw gas

- current status:
 - prototype at IFK is functional, but still requires frequent maintenance
 - design for low maintenance
TTA300 device available at Ratfish
 - next step: production of first serial device

FID Online Total Tar Analyzer TTA300



Analyzer measurement principle – loading phase



- Producer gas is sucked through sample loops via two venturi nozzles.
- Both loops filled successively.
- Detector signal is zero.

FID Online Total Tar Analyzer TTA300

Analyzer measurement principle – analyzing phase




- Volume of sample loops 1 & 2 flushed to FID successively.
- Tar substances of sample loop 1 are condensed on a cooled “tar filter”.
- Measurable components (groups):
 - Total hydrocarbon
 - Non-condensable HC
 - “Tars” (condensable HC)

FID Online Total Tar Analyzer TTA300

Comparative measurements with „Tar Protocol“

Gasification of wood pellets: 800°C, air ratio 0,3



Validation of measurement principle:

- Measurement cycle duration of 84 seconds corresponds to over 300 measurement cycles.
- Online measurements follow the wet chemical comparative measurements very clearly over total timeframe.
- Changes in tar concentrations by fluctuations in air ratio can be detected immediately.

ECN Classification:

2: Heterocyclic HC (e.g. pyridine, phenol, cresol)

3: 1-ring HC (e.g. xylene, styrene toluene)



4: 2, 3-ring PAH (e.g. naphthalene, anthracene, phenanthrene)

5: > 3-Ring PAH (e.g. fluoranthene, pyrene, chrysene)


<http://www.thersites.nl/classification.aspx>

FID Online Total Tar Analyzer TTA300

Comparative measurements with „Tar Protocol“



- Absolute deviations in comparative measurements also in other gasifier operating conditions in a range of max. $1\ g\ m^{-3}$.
- At *low* tar concentrations ($< 3\ g\ m^{-3}$) additional derivation is getting higher

 TTA 300 has good comparability to tar protocol and is an available device for raw gas monitoring


**sulfur measurement
at IFK**

tar protocol at IFK

Sewage sludge gasification 20 kW gasification tests

Tar and sulfur reduction with limestone as additive/sorbent


fuel: dried sewage sludge



analysis	in % kg/kg
H ₂ O	6.5
ash	47.6
C	51.0
H	6.9
O	32.0
N	7.5
S	2.4
Cl	0.2





Online gas measurement and sample gas preparation at IFK

Low cost H₂S and COS measurement via micro-GC



- gas washing with isopropanol and water necessary to protect online gas analyzers
- sour washing liquid enables H₂S and COS to pass the gas washing bottles
- H₂S and COS can be measured down to ~5 ppm


Sulfur reduction with limestone as sorbent




- higher limestone additive ratio reduces H₂S concentration
- H₂S concentration heads towards equilibrium for high additive ratios
- micro-GC als reliable sulfur measurement device

Tar protocol – wet chemical sampling

IFK's wet chemical tar sampling system




PSI's automatic liquid quench test at IFK's pilot plant




- Absorption liquid: Isopropanol
- Liquid is sampled and analysed
 - gravimetrically at IFK
 - GC-FID at FAU
 - GC-FID and GC-SCD at PSI → additional N and S-containing tars for residue gasification

Sar reduction with limestone as bed additive




tar classification




Class	anal.	description
ECN 1	gravimetric	gravimetric tars, condensable
ECN 2	GC-MS	water soluble tars
ECN 3	GC-MS	light aromatics, not condensable
ECN 4	GC-MS	light poly-aromatics (PAH), condensable
ECN 5	GC-MS	heavy PAHs, condensable

- substantial reduction of gravimetric tars
- higher limestone additive ratios needed for GC-MS tar reduction

- additional relevant components quantified by PSI with GC-FID/SCD:




Gravimetric tar composition for residue gasification



solid formation during gas sampling:

- ammoniumcarbonate:
 $2\text{NH}_3 + \text{CO}_2 + \text{H}_2\text{O} \rightarrow (\text{NH}_4)_2\text{CO}_3$
- ammoniumchloride:
 $\text{NH}_3 + \text{HCl} \rightarrow \text{NH}_4\text{Cl}$
- centrifugation removes most solids
- some are detected as grav. tar

elemental analysis of gravimetric tars



precipitation in tar sample

sample collection centrifugation → gravimetric analysis

**planned online GC for
gasification**

Online-GC for gasification

planned acquisition at IFK

requirements:

- sample line temperature incl. valves min. temp. 250°C
- measurement time < 30 min
- measured components
 - BTX (→ biggest fraction of „tars“)
 - C₁-C₆
 - NH₃
 - HCl
 - H₂O
 - permanent gases
- online coupling to IFK's gasification rigs → instant data

Summary and conclusion



- TTA 300 has good comparability to tar protocol and is an available device for raw gas monitoring
- H₂S and COS measurement with micro-GC is reliable for raw gas in residue gasification
- liquid sampling / tar protocol is still needed to characterize tars
 - for residue gasification additional GC detectable tar species need to be considered
 - solid formation needs to be considered for gravimetric analysis
 - endeavor to reduce liquid sampling effort by PSI's automatic liquid quench system
- new online GC planned at IFK



University of Stuttgart
Germany

Thank you!



M.Sc. Max Schmid

e-mail max.schmid@ifk.uni-stuttgart.de

phone +49 711 685-63394

fax +49 711 685-63491

University of Stuttgart
Institute of Combustion and Power Plant Technology
Pfaffenwaldring 23 • 70569 Stuttgart • Germany

